

To whom it may concern

RE: Aesculap C-Clamp for Unitrac Pneumatic Arm (RT079R) & Interface adapter (RT020R)

Please note that the Aesculap C-Clamp & interface adapter should be machine washed in a thermo washer disinfectant compatible to CFPP 01-01 Part D (HTM 2030) recommendations and should be subjected to a cleaning phase of 75 degrees Celsius and a disinfecting phase of between 90 to 95 degrees Celsius.

Chemistries used in this phase should be of a pH between 7.0 to 11.0. Please note, extended exposure to higher pH levels will result in damage to the surface of this device.

Sterilisation should follow CFPP 01-01 Part C (HTM 2010) recommendations and the C-Clamp & interface adapter can be sterilised between 134 to 137 degrees Celsius for a minimum holding time of three minutes to comply with the UK regulations

Aesculap develops and produces medical devices according to the European medical Device Directive (MDD 93/42/EEC). Our compliance with these regulations is demonstrated by displaying all of our devices with the CE mark, and all products and processes are strictly validated in accordance with MDD 93/44/EEC.

Amongst other things, MDD 93/44/EEC validates cleaning and reprocessing of medical devices, and Aesculap has proven effective reprocessing, with an internal process for every reusable medical device, including the C-Clamp & interface adapter for the Unitrac pneumatic arm.

If reprocessed in line with the above specifications and CFPP (HTM) recommendations, as validated, you can be assured that the C-Clamp & interface adapter will be fully decontaminated and sterile.

Please do not hesitate in contacting me directly should you require any additional information.

Kind Regards

Yours sincerely



Alan Barr
Clinical Manager Endosurgery
Aesculap Endoscopic Technology

Mobile : 07772 115856
email : allan.barr@bbraun.com