

## Be safe with SQ.line® – Probability of reprocessing rework is reduced!

Comparison of protein residues after washing – SQ.line vs. traditional design



SQ.line



Traditional instrument design



### Key Facts at a glance:

Comparing protein residues of double action instruments after Washer/Disinfector cycle, the **SQ.line design** offers advantages as follows:

- The **limit of protein residues** is reached at an **earlier stage** during cleaning and disinfection
- Therefore the risk of **reprocessing rework** is **reduced**



## Executive Summary - Test performance:

### Test purpose:

Prove that the new SQ.line double action instruments have a lower residual protein quantity compared to the traditional design for the corresponding cleaning time.

### Test object:


LX157NR vs. LX157R - defined as worst case design for this product group.

As competitors have equal constructions for double action instruments, the results are transferable to competitors products as well.

### Overview test procedure:

- 1<sup>st</sup> step: Basic cleaning of all instruments
- 2<sup>nd</sup> step: Soiling using sheep blood
- 3<sup>rd</sup> step: Cleaning
- 4<sup>th</sup> step: BCA method (a biochemical assay is used for quantitation of total protein in a sample) was used to determine the residual protein per test sample

### Soiling:



- After the complete soiling, the instrument was opened and closed 5 times.
- The instruments were dried for 1h at room temperature.

### Machine Cleaning:

- Test runs according to standard parameters of "Aesculap Validated Reprocessing Procedures AVA-V6" with interruption of the alkaline cleaning process after 5 and 3 minutes.
- After cleaning, the instruments were dried for 24h at room temperature.

### Cleaning Parameters:

| Step              | Temperature [°C] | Duration [min]   | Concentration [ml/l] | Medium  |
|-------------------|------------------|--|----------------------|---|
| Pre rinse         | 25               | 3  | -                    | deionised water                               |
| Alkaline cleaning | 55               | 1 <sup>st</sup> 5 minutes<br>2 <sup>nd</sup> 3 minutes | 5                    | Helimatic Cleaner Alkaline<br>deionised water |
| Final rinse       | cold             | 2  | -                    | deionised water                               |

**Note:** Alkaline cleaning time in daily CSSD routines takes 10 minutes. Objective of the above mentioned procedure is to proof, that SQ.line instruments reach the protein limits at an even earlier stage.

**Note:** The limit of protein residues for double action instruments is <100 µg.

### Test results:

|                           | Average absolute amount of protein (minus blank value) µg / sample 5 minutes* | Average absolute amount of protein (minus blank value) µg / sample 3 minutes* |
|---------------------------|---|---|
| LX157NR - SQ.line         | 37,87*  | 41,49*  |
| LX157R – traditional line | 63,93*  | 105,72*   |

\* Method-related error of ± 10µg / sample; two independent test runs

The SQ.line forceps show a 61% lower amount of protein residues after just 3 minutes of alkaline cleaning. Even after additional cleaning time of total 5 minutes the SQ.line forceps still show a 41% lower amount of protein residues than the traditional product line.

### Conclusion:

**A lower level of protein residues is reached at an earlier stage of the cleaning process compared to traditional instrument designs. Therefore the risk of reprocessing rework is reduced.**