

Omnifix® Luer Lock Solo (Type LDS)

3-Piece Fine Dosage Syringes with Luer Lock Connector



PRODUCT FEATURES

- Low Dead Space (LDS): Plunger with displacement spike **1**
- 0.01 / 0.1 ml graduation **2**
- Highly transparent barrel
- High contrast black graduation
- Integrated plunger backstop **3**
- Luer Lock connector
- Barrel made of polypropylene
- Plunger made of polystyrene
- Rubber piston with double sealing ring made of polyisoprene
- Not manufactured with PVC, Latex and DEHP
- Sterile single packaging
- Box of 100 pcs.

USER BENEFITS

- 0.01 / 0.1 ml graduation for precise dosage of smallest volumes
- Type LDS allows for minimal residual volume due to displacement spike
- Highly transparent barrel with high contrast black graduation for easy reading
- Rubber piston with double sealing ring allows slow aspiration and injection
- Safe plunger backstop to prevent unintentional plunger withdrawal

Omnifix® Luer Lock Solo (Type LDS), manufactured according to ISO 7886-1

Size	Cone	Graduation	Total Graduated Capacity	Units Per Box	Code No. (REF)
1 ml	Centric	0.01 ml	1.0 ml	100 pcs.	9167006V
2 ml	Centric	0.1 ml	2.5 ml	100 pcs.	4617029LDS