



# AESCULAP® PRIME*Line*® PRO

THE AESCULAP® STERILE CONTAINER SYSTEM

# THE AESCULAP® STERILE CONTAINER SYSTEM



The right INSTRUMENTS

At the right TIME

In sterile CONDITION

At the right PLACE OF USE

At a high level of PROFITABILITY

In a reliable PACKAGING SYSTEM




Fast, reliable and  
efficiently organized

## THE AESCULAP® STERILE CONTAINER SYSTEM

- Integrated permanent germ barrier system made of PTFE designed for 5,000 sterilisation cycles. No filter changes required - saves time and processing costs
- Stable, robust aluminium lid - reliable and long-lasting
- Outstanding cleaning results
- Good drying properties
- Available in five colours (red, blue, green, gold, silver) to help you organize your instrument sets





Article number	Size / colour	Dimensions mm (length x width x height)
JP101	1/1 Lid, red	588 x 281 x 36
JP102	1/1 Lid, blue	588 x 281 x 36
JP103	1/1 Lid, green	588 x 281 x 36
JP104	1/1 Lid, gold	588 x 281 x 36
JP105	1/1 Lid, silver	588 x 281 x 36
JP111	3/4 Lid, red	465 x 281 x 36
JP112	3/4 Lid, blue	465 x 281 x 36
JP113	3/4 Lid, green	465 x 281 x 36
JP114	3/4 Lid, gold	465 x 281 x 36
JP115	3/4 Lid, silver	465 x 281 x 36
JP121	1/2 Lid, red	298 x 281 x 36
JP122	1/2 Lid, blue	298 x 281 x 36
JP123	1/2 Lid, green	298 x 281 x 36
JP124	1/2 Lid, gold	298 x 281 x 36
JP125	1/2 Lid, silver	298 x 281 x 36
JP131	1/2 Lid for JK339, red	298 x 281 x 36
JP132	1/2 Lid for JK339, blue	298 x 281 x 36
JP133	1/2 Lid for JK339, green	298 x 281 x 36
JP134	1/2 Lid for JK339, gold	298 x 281 x 36
JP135	1/2 Lid for JK339, silver	298 x 281 x 36

Compatible with all bottom parts (1/1, 1/2, 3/4) manufactured since 1988  
 Consumables and accessories can be found in catalogue no. C40402



Fast, reliable and  
efficiently organized  
**THE AESCULAP®  
STERILE CONTAINER  
SYSTEM**

# AESCULAP® – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany  
Phone +49 7461 95-0 | Fax +49 7461 95-2600 | [www.aesculap.com](http://www.aesculap.com)

The main product trademark "Aesculap" and the product trademarks "PrimeLine" are registered trademarks of Aesculap AG. Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.